

Pipette Calibration Table of Z-Factors

Temperature (°C)	Z-Factor
15	1.0019
15.5	1.0020
16	1.0021
16.5	1.0022
17	1.0023
17.5	1.0024
18	1.0025
18.5	1.0026
19	1.0027
19.5	1.0028
20	1.0029
20.5	1.0030
21	1.0031
21.5	1.0032
22	1.0033
22.5	1.0034
23	1.0035
23.5	1.0036
24	1.0037
24.5	1.0038
25	1.0039

From: Pipette Performance Check. Wisconsin State Laboratory of Hygiene,
ESS INO GENOP 200, Revision 1, January, 2001.